

ELECTROSURGICAL  
INSTRUMENTS  
CATALOGUE



European Bipolar Forceps .....	01 - 07
European Irrigation Bipolar Forceps .....	08
European Non-Stick Bipolar Forceps .....	09 - 12
USA 2 Pin Bipolar Forceps .....	13 - 16
UAA 2-Pin Non-Stick Bipolar Forceps .....	17 - 19
Silver Glide Non-Stick Bipolar Forceps .....	20 - 22
Single Use USA 2-Pin Bipolar Forceps .....	23 - 25
Single Use European Bipolar Forceps .....	26 - 27
Monopolar Cables .....	28
Bipolar Forceps Cables .....	29
Diathermy Instruments .....	30 - 31
Bipolar Artery Sealer .....	32
Bipolar Electrodes .....	33
Electrodes 2.4mm .....	34
Conization Electrodes .....	35
Electrodes 4.0mm .....	36
Patient Return Plate .....	37
Laparoscopic Electrodes .....	38
Electrosurgical Instruments for Gynecology .....	39 - 42

McPherson (Str.)		8001
L= 9.0 cm		
Tip 0.5 mm		
McPherson (crv.)		8002
L= 11.5 cm		
Tip 0.5 mm		
Jeweler # 5		8003
L= 11.5 cm		
Tip 0.0 mm		
Adson		8003
L= 12.1 cm		
Tip 0.1 mm		
Jeweler # 4		8004
L= 11.5 cm		
Tip 0.5 mm		
Jeweler # 7		8005
L= 11.5 cm		
Tip 0.5 mm		
Coaptation		8006
L= 10.8 cm		
Tip 0.5 mm		

Coaptation Nadler		8007
L= 10.2 cm		
Tip 0.5 mm		
Iris (Str.)		8008
L= 10.2 cm		
Tip 0.5 mm		
Iris (crv.)		8009
L= 10.8 cm		
Tip 1.0 mm		
Straight		8010
L= 16.0 cm		
Tip 1.0 mm		
Straight		8011
L= 18.0 cm		
Tip 1.0 mm		
Straight		8012
L= 20.0 cm		
Tip 1.0 mm		
Straight		8013
L= 22.0 cm		
Tip 1.0 mm		

Straight (Ang.)		8014
L= 16.0 cm		
Tip 1.0 mm		

Straight (Ang.)		8015
L= 18.0 cm		
Tip 1.0 mm		

Straight (Ang.)		8016
L= 20.0 cm		
Tip 1.0 mm		

Straight (Ang.)		8017
L= 22.0 cm		
Tip 1.0 mm		

Bayonet		8018
L= 16.0 cm		
Tip 1.0 mm		

Bayonet		8019
L= 18.0 cm		
Tip 1.0 mm		

Bayonet		8020
L= 20.0 cm		
Tip 1.0 mm		

Bayonet (Ang.Up)

L= 16.0 cm

Tip 1.0 mm



8021

Bayonet (Ang.Up)

L= 18.0 cm

Tip 1.0 mm



8022

Bayonet (Ang.Up)

L= 20.0 cm

Tip 1.0 mm



8023

Bayonet (Ang.Down)

L= 16.0 cm

Tip 1.0 mm



8024

Bayonet (Ang.Up)

L= 18.0 cm

Tip 1.0 mm



8025

Bayonet (Ang.Up)

L= 20.0 cm

Tip 1.0 mm



8026

McPherson (Str.)

L= 9.0 cm

Tip 0.5 mm



8027

McPherson (crv.)		8028
L= 11.5 cm		
Tip 0.5 mm		
Jeweler # 5		8029
L= 11.5 cm		
Tip 0.5 mm		
Jeweler # 5		8030
L= 11.5 cm		
Tip 0.5 mm		
Jeweler # 7		8031
L= 11.5 cm		
Tip 0.5 mm		
Adson		8032
L= 12.1 cm		
Tip 1.0 mm		
Straight		8033
L= 16.0 cm		
Tip 1.0 mm		
Straight		8034
L= 18.0 cm		
Tip 1.0 mm		

Straight		8035
L= 20.0 cm		
Tip 1.0 mm		
Straight		8036
L= 22.0 cm		
Tip 1.0 mm		
Straight (Ang.)		8037
L= 16.0 cm		
Tip 1.0 mm		
Straight (Ang.)		8038
L= 18.0 cm		
Tip 1.0 mm		
Straight (Ang.)		8039
L= 20.0 cm		
Tip 1.0 mm		
Bayonet		8040
L= 16.0 cm		
Tip 1.0 mm		
Bayonet		8041
L= 18.0 cm		
Tip 1.0 mm		

Bayonet		8042
L= 22.0 cm		
Tip 1.0 mm		
Bayonet (Ang.Up)		8043
L= 16.0 cm		
Tip 1.0 mm		
Bayonet (Ang.Up)		8044
L= 18.0 cm		
Tip 1.0 mm		
Bayonet (Ang.Up)		8045
L= 20.0 cm		
Tip 1.0 mm		
Bayonet (Ang.Down)		8046
L= 16.0 cm		
Tip 1.0 mm		
Bayonet (Ang.Up)		8047
L= 18.0 cm		
Tip 1.0 mm		
Bayonet (Ang.Up)		8048
L= 20.0 cm		
Tip 1.0 mm		

Straight  
L= 16.0 cm  
Tip 1.0 mm



8049 IR

Straight  
L= 20.0 cm  
Tip 1.0 mm



8050 IR

Bayonet  
L= 16.0 cm  
Tip 1.0 mm



8051 IR

Bayonet  
L= 20.0 cm  
Tip 1.0 mm



8052 IR

McPherson (Str.)

L= 9.0 cm

Tip 0.5 mm



8053 NS

McPherson (crv.)

L= 11.5 cm

Tip 0.5 mm



8054 NS

Jeweler # 5

L= 11.5 cm

Tip 0.5 mm



8055 NS

Adson

L= 12.1 cm

Tip 1.0 mm



8056 NS

Jeweler # 4

L= 11.5 cm

Tip 0.5 mm



8057 NS

Jeweler # 7

L= 11.5 cm

Tip 0.5 mm



8058 NS

Coaptation

L= 10.8 cm

Tip 0.5 mm



8059 NS

Coaptation Nadler		8060 NS
L= 10.2 cm		
Tip 0.5 mm		
Iris (Str.)		8061 NS
L= 10.8 cm		
Tip 0.5 mm		
Iris (crv.)		8062 NS
L= 10.8 cm		
Tip 1.0 mm		
Straight		8063 NS
L= 18.0 cm		
Tip 1.0 mm		
Straight		8064 NS
L= 20.0 cm		
Tip 1.0 mm		
Straight		8065 NS
L= 20.0 cm		
Tip 1.0 mm		
Straight		8066 NS
L= 22.0 cm		
Tip 1.0 mm		

Straight (Ang.)

L= 16.0 cm

Tip 1.0 mm



8067 NS

Straight (Ang.)

L= 18.5 cm

Tip 1.0 mm



8068 NS

Straight (Ang.)

L= 20.0 cm

Tip 1.0 mm



8069 NS

Straight (Ang.)

L= 22.0 cm

Tip 1.0 mm



8070 NS

Bayonet

L= 16.0 cm

Tip 1.0 mm



8071 NS

Bayonet

L= 18.0 cm

Tip 1.0 mm



8072 NS

Bayonet

L= 20.0 cm

Tip 1.0 mm



8073 NS

Bayonet (Ang.Up)

L= 16.0 cm

Tip 1.0 mm



8074 NS

Bayonet (Ang.Up)

L= 18.0 cm

Tip 1.0 mm



8075 NS

Bayonet (Ang.Up)

L= 20.0 cm

Tip 1.0 mm



8076 NS

Bayonet (Ang.Down)

L= 16.0 cm

Tip 1.0 mm



8077 NS

Bayonet (Ang.Up)

L= 18.0 cm

Tip 1.0 mm



8078 NS

Bayonet (Ang.Up)

L= 20.0 cm

Tip 1.0 mm



8079 NS

McPherson (Str.)		8080
L= 9.0 cm		
Tip 0.5 mm		
McPherson (Crv.)		8081
L= 9.0 cm		
Tip 0.5 mm		
Coaptation Nadler		8082
L= 10.8 cm		
Tip 0.5 mm		
Coaptation		8083
L= 10.8 cm		
Tip 0.5 mm		
Jeweler # 5		8084
L= 11.0 cm		
Tip 0.5 mm		
Jeweler # 4		8085
L= 11.5 cm		
Tip 0.5 mm		
Tenzel		8087
L= 11.5 cm		
Tip 0.5 mm		

Jeweler # 7		8088
L= 11.5 cm		
Tip 0.5 mm		
Iris (Str.)		8089
L= 10.8 cm		
Tip 0.5 mm		
Iris (Crv.)		8090
L= 10.8 cm		
Tip 1.0 mm		
Adson		8091
L= 12.1 cm		
Tip 1.0 mm		
Semkin		8092
L= 14.0 cm		
Tip 1.0 mm		
Cushing		8093
L= 16.0 cm		
Tip 1.0 mm		
Cushing		8094
L= 11.5 cm		
Tip 1.0 mm		

Cushing (Gursh)		8095
L= 18.0 cm		
Tip 1.0 mm		
Gerald (Str)		8096
L= 17.8 cm		
Tip 1.0 mm		
Gerald (Crv)		8097
L= 17.8 cm		
Tip 1.0 mm		
Gerald Bayonet		8098
L= 17.8 cm		
Tip 1.0 mm		
Jansen Bayonet		8099
L= 15.0 cm		
Tip 1.0 mm		
Jansen Bayonet		8100
L= 18.0 cm		
Tip 1.0 mm		
Malis		8101
L= 17.8 cm		
Tip 1.0 mm		

Scollive Green-Wood		8102
L= 19.1 cm		
Tip 1.5 mm		
Yasargil		8103
L= 19.8 cm		
Tip 1.0 mm		
Hardy Bayonet		8104
L= 12.5 cm		
Tip 1.0 mm		
Hardy Bayonet (Ang)		8105
L= 21.5 cm		
Tip 1.0 mm		

McPherson (Str)

L= 9.0 cm

Tip 0.5 mm



8106 NS

McPherson (Crv)

L= 9.0 cm

Tip 0.5 mm



8107 NS

Jewelery #5

L= 11.5 cm

Tip 0.5 mm



8108 NS

Jeweler # 4

L= 11.5 cm

Tip 0.5 mm



8109 NS

Adson

L= 12.1 cm

Tip 1.0 mm



8110 NS

Semkin

L= 14.0 cm

Tip 1.0 mm



8111 NS

Cushing

L= 16.0 cm

Tip 1.0 mm



8112 NS

Cushing		8113 NS
L= 18.0 cm		
Tip 1.0 mm		
Gerald Bayonet		8114 NS
L= 17.8 cm		
Tip 1.0 mm		
Jansen Bayonet		8115 NS
L= 18.0 cm		
Tip 1.0 mm		
Scollive Green-Wood		8116 NS
L= 19.1 cm		
Tip 1.5 mm		
Yasargil		8117 NS
L= 19.8 cm		
Tip 1.0 mm		
Hardy Bayonet		8118 NS
L= 21.5 cm		
Tip 1.0 mm		

Semkin Irrigation		8119 IR
Gerald Bayonet Irrigation		8120 IR
Jansen Bayonet Irrigation		8121 IR
Scollive Green-Wood Irrigation		8122 IR
Hardy Bayonet Irrigation		8123 IR

McPherson (Str.)		8124 FNS
L= 9.0 cm		
Tip 0.5 mm		
McPherson (Crv.)		8125 FNS
L= 9.0 cm		
Tip 0.5 mm		
Jeweler # 5		8126 FNS
L= 11.5 cm		
Tip 0.5 mm		
Adson		8127 FNS
L= 12.1 cm		
Tip 1.0 mm		
Cushing (Str.)		8128 FNS
L= 16.0 cm		
Tip 1.0 mm		
Jansen Bayonet		8129 FNS
L= 18.0 cm		
Tip 1.0 mm		
Jansen Bayonet		8130 FNS
L= 20.0 cm		
Tip 1.0 mm		
Bayonet Full Non-Stick		8131 FNS
L= 16.0 cm		
Tip 1.0 mm		

Bayonet Full Non-Stick

L= 18.0 cm

Tip 1.0 mm



8132 FNS

Bayonet Full Non-Stick

L= 20.0 cm

Tip 1.0 mm



8133 FNS

Bayonet Full Non-Stick

L= 22.0 cm

Tip 1.0 mm



8134 FNS

Bayonet (Irrigation)

L= 20.0 cm

Tip 1.0 mm



8135 FNS

Bayonet (Irrigation)

L= 22.0 cm

Tip 1.0 mm



8136 FNS

Bayonet Full Non-Stick

L= 22.0 cm

Tip 1.0 mm



8137 FNS

Bayonet Full Non-Stick

L= 22.0 cm

Tip 0.5 mm



8138 FNS

Bayonet Full Non-Stick

L= 24.0 cm

Tip 1.0 mm



8139 FNS

Jeweler # 5 L= 11.5 cm Tip 0.5 mm		8140 FNS
Jeweler # 5 L= 11.5 cm Tip 0.5 mm		8142 FNS
Straight Full Non-Stick L= 12.0 cm Tip 1.0 mm		8143 FNS
Straight (Ang.) L= 16.0 cm Tip 1.0 mm		8144 FNS
Bayonet Full Non-Stick L= 18.0 cm Tip 1.0 mm		8145 FNS
Bayonet Full Non-Stick L= 20.0 cm Tip 1.5 mm		8146 FNS
Bayonet Full Non-Stick L= 22.0 cm Tip 1.5 mm		8147 FNS
Bayonet Full Non-Stick L= 22.0 cm Tip 1.0 mm		8148 FNS

Gerald Straight  
Single Use



8149 S

Jansen Bayonet  
Irrigation



8150 IRS

Jansen Bayonet  
Single Use



8151 S

Jeweler # 5  
Single Use



8152 S

Adson Forceps  
Single Use



8153 S

Cushing Forceps  
Size: 16.0 cm  
Tip: 1.0 mm



8154 S

Jansen Bayonet ang-up  
Size: 20.0 cm  
Tip: 1.0 mm



8154 S

Jeweler # 5 L= 11.5 cm Tip 0.5 mm		8155 FNS
Jeweler # 5 L= 11.5 cm Tip 0.5 mm		8156 FNS
Straight Full Non-Stick L= 12.0 cm Tip 1.0 mm		8157 FNS
Straight (Ang.) L= 16.0 cm Tip 1.0 mm		8158 FNS
Bayonet Full Non-Stick L= 18.0 cm Tip 1.0 mm		8159 FNS
Bayonet Full Non-Stick L= 20.0 cm Tip 1.5 mm		8160 FNS
Bayonet Full Non-Stick L= 22.0 cm Tip 1.5 mm		8161 FNS
Bayonet Full Non-Stick L= 22.0 cm Tip 1.0 mm		8162 FNS

Jansen Bayonet  
Size: 18.0 cm  
Tip: 1.0 mm



8163 S

Jeweler # 4  
Size: 11.5 cm  
Tip: 0.5 mm



8164 S

Adson Forceps  
Size: 12.1 cm  
Tip: 1.0 mm



8165 S

Jansen Bayonet  
Size: 18.0 cm  
Tip: 1.0 mm



8166 S

Jansen Bayonet  
Size: 20.0 cm  
Tip: 1.0 mm



8167 S

Bayonet Forceps

Size: 18.0 cm

Tip: 1.0 mm



8173 S

Bayonet Forceps

Size: 20.0 cm

Tip: 1.0 mm



8174 S

Straging Forceps

Size: 20.0 cm

Tip: 1.0 mm



8175 S

Straight Forceps (ang)

Size: 14.0 cm

Tip: 1.0 mm



8176 S

Straight Forceps

Size: 14.0 cm

Tip: 1.0 mm



8177 S

Jeweler # 5  
Single Use



8168 S

Straight Forceps  
Single Use



8169S

Bayonet Forceps  
Single Use



8170 S

Straight Forceps  
Single Use



8171 S

Straight Forceps (ang)  
Single Use



8172 S



**Silicone Monopolar Cables 10" (3.0m)**  
 8185  
 Instrument Connector: Screw Type  
 Generator Connector: Standard Bovie 8.0mm Plug  
 Compatibility: Valleylab, Conmed, Bovie



**Silicone Monopolar Cables 10" (3.0m)**  
 8186  
 Instrument Connector: 4.0mm Male Fitting  
 Generator Connector: Standard 4.0mm Banana Plug  
 Compatibility: Valleylab, Conmed & All ES Generators



**Silicone Monopolar Endoscopic Cable 10" (3.0m)**  
 8187  
 Instrument Connector: 4.0mm Female Plug  
 Generator Connector: Standard Bovie 8.0mm Plug  
 Compatibility: Valleylab, Conmed, Bovie



**Silicone Monopolar Endoscopic Cable 10" (3.0m)**  
 8188  
 Instrument Connector: 4.0mm Female Plug  
 Generator Connector: Standard 4.0mm Banana Plug  
 Compatibility: Valleylab, Conmed & All ES Generators



**Silicone Diathermy Cable 10" (3.0m)**  
 8189  
 Instrument Connector: 4.8mm Male Fitting  
 Generator Connector: Standard Bovie 8.0mm Plug  
 Compatibility: Valleylab, Conmed, Bovie



**Silicone Diathermy Cable 10" (3.0m)**  
 8190  
 Instrument Connector: 4.8mm Male Fitting  
 Generator Connector: Standard 4.0mm Banana Plug  
 Compatibility: Valleylab, Conmed & All ES Generators



**Silicone Bipolar Cables 10" (3.0m)**  
 8178  
 Instrument Connector: USA 2-Pin Plug  
 Generator Connector: 4.0mm Banana Plug  
 Compatibility: Valleylab, Conmed, Codman, Olympus and Mbst Bipolar Generator



**Silicone Bipolar Cables 10" (3.0m)**  
 8179  
 Instrument Connector: European Flat Plug  
 Generator Connector: Martin Type  
 Compatibility: Aesculap, Martin, Wolf, Berchold



**Silicone Bipolar Cables 10" (3.0m)**  
 8180  
 Instrument Connector: European Flat Plug  
 Generator Connector: 4.0mm Banana Plug  
 Compatibility: Valleylab, Conmed, Codman, Olympus & Mbst Bipolar Generator



**Silicone Bipolar Cables 10" (3.0m)**  
 8181  
 Instrument Connector: USA2-Pin Plug  
 Generator Connector: Martin Type  
 Compatibility: Aesculap, Martin, Wolf, Berchold



**Silicone Bipolar Cables 10" (3.0m)**  
 8182  
 Instrument Connector: USA2-Pin Plug  
 Generator Connector: Erbe Type  
 Compatibility: Eibe, Siemens, Downs



**Silicone Bipolar Cables 10" (3.0m)**  
 8183  
 Instrument Connector: USA2-Pin Plug  
 Generator Connector: Pin Distance 29.0mm  
 Compatibility: Valleylab, EMS, Dolley

**Silicone Bipolar Cables 10" (3.0m)**  
 8184  
 Instrument Connector: European Flat Plug  
 Generator Connector: Pin Distance 29.0mm  
 Compatibility: Valleylab, EMS, Dolley





**Turnur Warwick Scissor**  
(Straight)  
8191



**Turnur Warwick Scissor**  
(Curved)  
8192



**Surgical Scissor**  
(Straight)  
8193



**Surgical Scissor**  
(Curved)  
8194



**Abbey Sub Needle**  
17.5cm  
8195



**Willson Hey Forceps**  
(Straight)  
8196

**Richie Forceps**  
(Straight)  
8197



**Willson Hey Forceps**  
(Curved)  
8198



**Artery Forceps**  
(Straight)  
8199



**Artery Forceps**  
(Curved)  
8200



Adson		8201
Mcdone		8202
Gillie		8203
Dissecting Forceps		8204
Gerald (1x2)		8205
Adson Bayonet		8206
Adson Bayonet Angle Up		8207
Remington Hob		8208



Artery Forceps Bi-Clamp

8209



Artery Forceps Non Stick

8210




Artery Forceps

8211



Artery Forceps

8212



Silicone Bi-Clamp Connecting cord



**Straight**  
8213

**Angled**  
8214

**Angled**  
8215

**Angled 45°**  
8216

**Angled 90°**  
8217



**Bi-Turbinate Probe**  
Size: 17.5cm  
8218



**Blade Electrode**  
Disposable: 8219



**Ball Electrode (crv.) (3.0 mm)**  
Disposable: 8222



**Needle Electrode**  
Disposable: 8225



**Needle Electrode (crv.)**  
Disposable: 8228



**Ball Electrode (2.0 mm)**  
Disposable: 8231




**Ball Electrode (crv.) (2.0 mm)**  
Disposable: 8233



**Bladeelectrode (Shaft 13 .0 cm )**  
Disposable: 8236



**Bladeelectrode (Shaft 13 .0 cm )**  
Disposable: 8238



**Bladeelectrode (Shaft 13 .0 cm )**  
Disposable: 8240



**Ball Electrode (3.0 mm)**  
Disposable: 8220



**Blade Electrode (crv.)**  
Disposable: 8223



**Ball Electrode (4.0 mm)**  
Disposable: 8226



**Ball Electrode (crv.) (4.0 mm)**  
Disposable: 8230



**Ball Electrode (5.0 mm)**  
Disposable: 8232



**Ball Electrode (crv.) (5.0 mm)**  
Disposable: 8234



**Circle Loop Electrode 4.0 mm**  
Reusable: 8221



**Circle Loop Electrode 4.0 mm**  
Reusable: 8224



**Circle Loop Electrode 8.0 mm**  
Reusable: 8227



**Diamond Electrode 8.0 mm**  
Reusable: 8229



**Tip Cleaner** 8235



**Holster** 8237



**Bladeelectrode (Shaft 13 .0 cm )**  
Disposable: 8239




**Conizationelectrode**  
Diposable : 8241



Conization Electrode  
Disposable 8242



Conization Electrode  
Disposable 8243



Conization Electrode  
Disposable 8244



Conization Electrode  
Disposable 8245



Conization Electrode  
Disposable 8246



Extender Electrode  
8247 Size: 7.0 cm



Extender Electrode  
8248 Size: 13.0 cm



Extender Electrode 8249 Size: 18.0 cm



**Blade Electrode**  
Reusable 8250



**Blade Electrode (Crv)**  
Reusable 8252



**Needle Electrode**  
Reusable 8254



**Needle Electrode (Crv)**  
Reusable 8257



**Ball Electrode (2mm)**  
Reusable 8259



**Ball Electrode (Crv) (2mm)**  
Reusable 8262



**Blade Electrode (Shaft 13.0cm)**  
Reusable SM-8265



**Needle Electrode (Shaft 13.0cm)**  
Reusable 8267



**Ball Electrode 3mm (Shaft 13.0 cm)**  
Reusable 8268



**Ball Electrode 4mm (Shaft 13.0cm)**  
Reusable 8270



**Ball Electrode 5mm (Shaft 13.0cm)**  
Reusable 8271



**Ball Electrode**  
Reusable 8251



**Ball Electrode (Crv) 3mm**  
Reusable 8253



**Ball Electrode (4mm)**  
Reusable 8255



**Ball Electrode (Crv) (4mm)**  
Reusable 8258



**Ball Electrode (5mm)**  
Reusable 8260



**Ball Electrode (Crv) (5mm)**  
Reusable 8263



**Round Loop Electrode (7.0cm) 4.0mm**  
Reusable 8256



**Round Loop Electrode (7.0cm) 6.0mm**  
Reusable 8259



**Round Loop Electrode (7.0cm) 8.0mm**  
Reusable 8261



**Round Loop Electrode (7.0cm) 10.0mm**  
Reusable 8264



**Plate Electrode Small**  
Reusable 8266



**Plate Electrode Large**  
Reusable 8269

**Patient Return Plate Reusable**  
5m Cable 8271



**Patient Return Plate Reusable**  
(5m Cable) 8272



**Patient Return Cable**  
8273- a(5m cable)  
8274- b(3m Cable)



**Patient Return Cable**  
8275- a(5m cable)  
8276- b(3m Cable)



**S.S Patient Return Plate**  
(16x12cm) 8277

# Laparoscopic Electrodes

## WITHOUT SUCTION



8278-Hook Electrode (36cm)



8280-Hook Electrode (36cm)



8282-Needle-Hook Electrode (36cm)



8284-Spatula-Hook Electrode (36cm)



8286-Ball-Hook Electrode (36cm)

## WITH SUCTION



8279L-Hook Electrode (Irrigation) (36cm)



8281J-Hook Electrode (Irrigation) (36cm)



8283-Needle-Hook Electrode (Irrigation) (36cm)



8285-Spatula-Hook Electrode (Irrigation) (36cm)



8287-Ball-Hook Electrode (Irrigation) (36cm)

## Laparoscopic Electrodes





**Graves Vaginal Speculum**  
with smoke evacuation tube  
**Small:** 8288  
**Medium:** 8289  
**Large:** 8290



**VU-Max Speculum**  
with smoke evacuation tube  
**Small:** 8291  
**Medium:** 8292  
**Large:** 8293



**Pederson Vaginal Speculum**  
with smoke evacuation tube  
**Small:** 8294  
**Medium:** 8295  
**Large:** 8296



**VU-More (Graves) Speculum**  
with smoke evacuation tube  
**Small:** 8300  
**Medium:** 8301  
**Large:** 8302

**Cusco Speculum**  
with smoke evacuation tube  
**Small:** 8297  
**Medium:** 8298  
**Large:** 8299



**Collin Speculum**  
with smoke evacuation tube  
**Small:** 8303  
**Medium:** 8304  
**Large:** 8305





**Weisman Speculum with Suction (Left Open)**  
**Small:** 8306  
**Medium:** 8307  
**Large:** 8308



**Weisman Speculum with Suction (Right Open)**  
**Small:** 8309  
**Medium:** 8310  
**Large:** 8311



**Dressing Forceps**  
8312(21.0 cm)  
8313(25.0 cm)



**Tissue Forceps**  
8314(21.0 cm)  
8315 (25.0 cm)



**Iris Hook**  
8316(25.0 cm)



**Emmet Hook**  
8317(25.0 cm)



**Straight Hook**  
8318(25.0 cm)



**Two-Prong Hook**  
8319(25.0 cm)





Lateral Vaginal  
Retractor 8320



Kogan Endocervical  
Speculum 8321



Lateral Vaginal Retractor  
(wide-Open) 8322



Rigby Retractor  
8323




Rigby Retractor  
8324





**Kogan Endocervical Speculum**  
8329



**Lateral Vaginal Retractor**  
8328



**Graves Speculum with smoking evacuation tube**  
Small: 8330  
Medium: 8331  
Large: 8332



**Pederson Graves Speculum with smoking evacuation tube**  
Small: 8336  
Medium: 8337  
Large: 8338

**Collin Speculum with smoking evacuation tube**  
Small: 8333  
Medium: 8334  
Large: 8335





1800 W HAWTHORNE LN STE P2 WEST  
CHICAGO, IL 60185  
Phone: (877) 400-2601  
Email: [info@ameridentinc.com](mailto:info@ameridentinc.com)  
Web: [www.ameridentinc.com](http://www.ameridentinc.com)